

## Anti-Syk antibody (Central Region) [M373] (STJA0041385)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

### PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A Purified
Dilution Range	WB 1:250
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN <sub>3</sub> and 50% glycerol
Isotype	IgG1
Molecular Weight	72
Storage Instruction	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.

### TARGET INFORMATION

Gene ID	<a href="#">6850</a>
Gene Symbol	<a href="#">SYK</a>
UniProt ID	<a href="#">KSYK_HUMAN</a>
Immunogen Region	Central Region
Immunogen Sequence	Human
Specificity	This antibody detects a 72 kDa* protein corresponding to the molecular mass of Syk on SDS-PAGE immunoblots of adult mouse spleen.

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



anti-Syk antibody Western blot analysis adult mouse spleen lysate STJA0041385