

Anti-Phospho-c-Src-Tyr-215 antibody (STJA0041394)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:1000
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG
Molecular Weight	60
Storage Instruction	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.

TARGET INFORMATION

Gene ID	6714
Gene Symbol	SRC
UniProt ID	SRC_HUMAN
Immunogen Sequence	Human
Specificity	This antibody was cross-adsorbed to phospho-tyrosine coupled to agarose then affinity purified using phospho-c-Src (Tyr-215) peptide (without carrier). The antibody detects a 60 kDa* protein on SDS-PAGE immunoblots of mouse SYF cells transformed with

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot c-Src phosphorylation in mouse SYF cells with pervanadate treatment STJA0041394