

Anti-Semaphorin-3E antibody (N-Term) (STJA0041405)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB
Host / Source	Rabbit
Reactivity	Chicken/Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:1000
Formulation	PBS + 0.05% NaN ₃
Isotype	IgG
Molecular Weight	61/87
Storage Instruction	Recommended that the undiluted antibody be aliquoted into smaller working volumes (10-30 uL/vial depending on usage) upon arrival and stored long term at -20°C C or -80°C C, while keeping a working aliquot stored at 4°C C for short term. Avoid freeze/t

TARGET INFORMATION

Gene ID	9723
Gene Symbol	SEMA3E
UniProt ID	SEM3E_HUMAN
Immunogen Region	N-Term
Immunogen Sequence	Human
Specificity	The antibody detects 87 and 61 kDa bands corresponding to the full-length and furin-processed forms of Sema-3E in human PC-3 and rat PC12 cells.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot rat PC12 and human PC-3 cells using anti-Sema-3E antibody with blocking peptide STJA0041405