

Anti-Beta III-Tubulin antibody (C-Term) (STJA0041422)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IHC
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:2000 IHC 1:500 ICC 1:200
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN3 and 50% glycerol
Isotype	IgG
Molecular Weight	50
Storage Instruction	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.

TARGET INFORMATION

Gene ID	10381
Gene Symbol	TUBB3
UniProt ID	TBB3_HUMAN
Immunogen Region	C-Term
Immunogen Sequence	Human
Specificity	The antibody detects a 50 kDa* protein corresponding to the molecular mass of Beta III-Tubulin on SDS-PAGE immunoblots of purified brain tubulin and mouse brain tissue.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---