

Anti-Beta-Tubulin antibody (Central Region) (STJA0041424)

GENERAL INFORMATION

| | |
|---------------|-------------------------|
| Product Type | Primary antibodies |
| Applications | ELISA/WB/ICC |
| Host / Source | Rabbit |
| Reactivity | Chicken/Human/Mouse/Rat |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Clonality | Polyclonal |
| Conjugation | Unconjugated |
| Purification | Antigen Affinity Purified |
| Dilution Range | WB 1:1000 ICC 1:50 |
| Formulation | PBS + 1 mg/ml BSA, 0.05% NaN3 and 50% glycerol |
| Isotype | IgG |
| Molecular Weight | 50 |
| Storage Instruction | Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C. |

TARGET INFORMATION

| | |
|--------------------|---|
| Gene ID | 10381 |
| Gene Symbol | TUBB3 |
| UniProt ID | TBB3_HUMAN |
| Immunogen Region | Central Region |
| Immunogen Sequence | Human |
| Specificity | This antibody detects a 50 kDa* protein corresponding to the molecular mass of Beta-Tubulin on SDS-PAGE immunoblots of purified brain tubulin and mouse brain tissue. |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|