

Anti-Phospho-Alpha 6-Tubulin-Ser-165 antibody (STJA0041427)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB
Host / Source	Rabbit
Reactivity	Chicken/Human/Mouse/Rat/Xenopus

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:1000
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN3 and 50% glycerol
Isotype	IgG
Molecular Weight	50
Storage Instruction	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.

TARGET INFORMATION

Gene ID	84790
Gene Symbol	TUBA1C
UniProt ID	TBA1C_HUMAN
Immunogen Sequence	Human
Specificity	This antibody was cross-adsorbed to unphosphorylated Alpha 6-Tubulin (a.a. 160-169) peptide before affinity purification using Alpha 6-Tubulin (Ser-165) peptide (without carrier). The antibody detects a 50 kDa* protein corresponding to the molecular

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---