

Anti-Vav2 antibody (STJA0041438)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IP
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:1000 ICC 1:50 IP 1:100
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG
Molecular Weight	100
Storage Instruction	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.

TARGET INFORMATION


Gene ID	22325
Gene Symbol	Vav2
UniProt ID	VAV2_MOUSE
Immunogen Sequence	Mouse
Specificity	The antibody detects a 100 kDa* band corresponding to Vav2 on SDS-PAGE immunoblots of human A431 cells and detects several Vav isoforms in Jurkat cells.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot EGF-treated A431 cells probed with anti-Vav and anti-Vav2 antibodies STJA0041438



ICC staining Vav protein in rabbit spleen fibroblasts using antibody with Cy3 detection STJA0041438