

Anti-Vav3 antibody (STJA0041442)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:1000
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG
Molecular Weight	98
Storage Instruction	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.

TARGET INFORMATION

Gene ID	10451
Gene Symbol	VAV3
UniProt ID	VAV3_HUMAN
Immunogen Sequence	Human
Specificity	The antibody detects a 98 kDa* protein corresponding to the molecular mass of Vav3 on SDS-PAGE immunoblots of human Jurkat cells. The antibody also strongly detects an unidentified protein at 120 kDa in Jurkat and K562.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot Vav2 and Vav3 detection in JurkatHUVECand A431 cell lysates STJA0041442