

Anti-WAVE2 antibody (Central Region) (STJA0041449)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:1000 ICC 1:250
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG
Molecular Weight	84
Storage Instruction	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.

TARGET INFORMATION


Gene ID	10163
Gene Symbol	WASF2
UniProt ID	WASF2_HUMAN
Immunogen Region	Central Region
Immunogen Sequence	Human
Specificity	The antibody detects a 84 kDa* protein corresponding to the molecular mass of WAVE2 on SDS-PAGE immunoblots of human K562, mouse C2C12, and rat PC12 cells.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



anti-WAVE2 antibody WB analysis rat PC12 lysate with blocking peptide control STJA0041449



ICC staining WAVE2 in rat PC12 cells using antibody with Cy3 detection STJA0041449