

Anti-N-WASP antibody (N-Term) (STJA0041451)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:1000
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG
Molecular Weight	65
Storage Instruction	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.

TARGET INFORMATION

Gene ID	8976
Gene Symbol	WASL
UniProt ID	WASL_HUMAN
Immunogen Region	N-Term
Immunogen Sequence	Human
Specificity	This antibody detects a 65 kDa* protein corresponding to the molecular mass of N-WASP on SDS-PAGE immunoblots of neonatal rat brain lysate. It is also detects 65 kDa* proteins in A431, human endothelial, and SKN-SH cells. It does not recognize the 63

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot N-WASP phosphorylated and unphosphorylated forms in rat brain lysate STJA0041451