

Anti-Unphospho-N-WASP antibody (STJA0041454)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:1000 ICC 1:50
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG
Molecular Weight	65
Storage Instruction	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.

TARGET INFORMATION


Gene ID	8976
Gene Symbol	WASL
UniProt ID	WASL_HUMAN
Immunogen Sequence	Human
Specificity	This antibody was cross-adsorbed to phospho-N-WASP (S484/S485) peptide (without carrier) then affinity purified using unphosphorylated N-WASP (S484/S485) peptide (without carrier). The antibody detects a 65 kDa* protein on SDS-PAGE immunoblots of rat

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot anti-N-WASP antibody detection in neonatal rat brain lysate STJA0041454



ICC staining N-WASP phosphorylation in rabbit spleen fibroblasts with Cy3 detection STJA0041454