

## Anti-P2Y Purinoceptor 2 antibody (N-Term) (STJA0041458)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

### PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:500-1:1000 IHC 1:100-1:200
Formulation	PBS + 0.02% NaN <sub>3</sub>
Isotype	IgG
Storage Instruction	Recommended that the undiluted antibody be aliquoted into smaller working volumes (10-30 uL/vial depending on usage) upon arrival and stored long term at -20°C C or -80°C C, while keeping a working aliquot stored at 4°C C for short term. Avoid freeze/t

### TARGET INFORMATION

Gene ID	<a href="#">5029</a>
Gene Symbol	<a href="#">P2RY2</a>
UniProt ID	<a href="#">P2RY2_HUMAN</a>
Immunogen Region	N-Term
Immunogen Sequence	Human
Specificity	Reacts with the N-terminal region of human, mouse and rat P2Y2 protein.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---