

Anti-CNP antibody (Full-Length) (STJA0041459)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/IHC
Host / Source	Goat
Reactivity	Bovine/Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 3.513888889 IHC 2.125
Formulation	PBS + 50% glycerol and 5 mM NaN3
Isotype	IgG
Molecular Weight	46
Storage Instruction	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.

TARGET INFORMATION

Gene ID	1267
Gene Symbol	CNP
UniProt ID	CN37_HUMAN
Immunogen Region	Full-Length
Immunogen Sequence	Human
Specificity	Specific for endogenous levels of the ~46 kDa CNP protein.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---