

Anti-Adiponectin/Acrp30 antibody (15-244aa) [1G12] (STJA0041460)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | ELISA/WB |
| Host / Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES


| | |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------|
| Clonality | Monoclonal |
| Concentration | 1 mg/mL |
| Conjugation | Unconjugated |
| Purification | By protein-G affinity chromatography |
| Formulation | Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol |
| Isotype | IgG1k |
| Storage Instruction | For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles. |

TARGET INFORMATION

| | |
|------------------|-----------------------------|
| Gene ID | 9370 |
| Gene Symbol | ADIPOQ |
| UniProt ID | ADIPO_HUMAN |
| Immunogen Region | 15-244aa |

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot of human adiponectin Acrp30 antibody detecting full-length and globular domain