

Anti-Adiponectin/Acrp30 antibody (15-244aa) [5H7] (STJA0041461)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-G affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG1k
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	9370
Gene Symbol	ADIPOQ
UniProt ID	ADIPO_HUMAN
Immunogen Region	15-244aa

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC IF of Adiponectin in Hep3B cells stained with STJA0041461 antibody Alexa Fluor 488 secondary DAPI nucleus



Western blot of human Adiponectin Acrp30 antibody STJA0041461 showing 293T cell lysates with Adiponectin domains