

Anti-5'-Nucleotidase/CD73 antibody (27-252aa) [AT89E9] (STJA0041474)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF/FACS
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-G affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG2bk
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	4907
Gene Symbol	NT5E
UniProt ID	5NTD_HUMAN
Immunogen Region	27-252aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC IF of 5'-Nucleotidase CD73 in HeLa cells stained with STJA0041474 antibody



Western blot of human 5'-Nucleotidase CD73 antibody STJA0041474 showing protein expression in HepG2 WtDr and A431 cell lysates



Flow cytometry of CD73 5'-Nucleotidase in A431 cells with STJA0041474 antibody