

Anti-FABP1/L-FABP antibody (1-127aa) [2G4] (STJA0041489)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF/IHC/FACS
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-G affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG1k
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	2168
Gene Symbol	FABP1
UniProt ID	FABPL_HUMAN
Immunogen Region	1-127aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




ICC IF of FABP1 in Hep3B cells stained with STJA0041489 antibody Alexa Fluor secondary green and DAPI blue



Western blot of human FABP1 L-FABP antibody showing recombinant protein and cell lysates



Flow cytometry of FABP1 in Hep3B cells using Human FABP1 L-FABP antibody with Alexa Fluor 488 detection



Human FABP1 L-FABP antibody IHC staining in paraffin-embedded colon cancer tissue