

Anti-GSK-3 beta antibody (341-420aa) [1F7] (STJA0041493)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF/FACS
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-G affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG1k
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	2932
Gene Symbol	GSK3B
UniProt ID	GSK3B_HUMAN
Immunogen Region	341-420aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC IF of GSK-3 beta in HeLa cells with Alexa Fluor 488 secondary antibody and DAPI nuclear staining



Western blot of human GSK-3 beta protein in HeLa cell lysate using anti-human GSK3 beta antibody



Flow cytometry of GSK-3 beta in HeLa cells using STJA0041493 antibody with Alexa Fluor 488 detection