

Anti-HSP60 antibody (1-573aa) [2E4] (STJA0041505)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF/IHC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-G affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG1k
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3329
Gene Symbol	HSPD1
UniProt ID	CH60_HUMAN
Immunogen Region	1-573aa

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC IF of HSP60 in HeLa cells stained with STJA0041505 antibody Alexa Fluor 488 secondary DAPI nuclear counterstain



Western blot of Human HSP60 antibody in 293T HeLa and A549 cell lysates



Human HSP60 antibody IHC staining in paraffin-embedded colon cancer tissue with DAB detection