

Anti-MIF antibody (1-115 aa) [4E4] (STJA0041518)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-G affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG1k
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	4282
Gene Symbol	MIF
UniProt ID	MIF_HUMAN
Immunogen Region	1-115 aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




ICC IF of MIF in Balb 3T3 cells stained with STJA0041518 antibody



ICC IF of MIF in HeLa cells stained with human MIF antibody STJA0041518 and Alexa Fluor 488



Western blot of human MIF antibody STJA0041518 in Jurkat THP-1 HeLa and U937 cell lysates



Western blot showing human MIF antibody detection in HL-60 cell lysates using anti-mouse HRP secondary