

Anti-Nanog antibody (1-154aa) [5A10] (STJA0041519)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/IHC/FACS
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-G affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG2ak
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	79923
Gene Symbol	NANOG
UniProt ID	NANOG_HUMAN
Immunogen Region	1-154aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot of human Nanog antibody STJA0041519 in NIH3T3 cell lysate



Flow cytometry of Nanog in Hep3B cells stained with human Nanog antibody and Alexa Fluor 488



Human Nanog antibody IHC staining in paraffin-embedded seminoma tissue with DAB detection