

Anti-PP2C gamma/PPM1G antibody (317-546aa) [k1G6] (STJA0041534)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF/IHC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-G affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG1k
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	5496
Gene Symbol	PPM1G
UniProt ID	PPM1G_HUMAN
Immunogen Region	317-546aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC IF of PPM1G in 293T cells with STJA0041534 antibody Alexa Fluor 488 green and DAPI blue



Western blot of human PPM1G protein in 293T cell lysate using anti-PPM1G antibody



IHC staining of human breast cancer tissue with anti-PPM1G antibody using DAB detection