

Anti-RASSF1A antibody (1-340aa) [3F3] (STJA0041537)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-A affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG1k
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	11186
Gene Symbol	RASSF1
UniProt ID	RASF1_HUMAN
Immunogen Region	1-340aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




RASSF1A immunofluorescence in HeLa cells with Alexa Fluor 488 and DAPI staining



ICC IF of RASSF1A in A431 cells with STJA0041537 antibody and Alexa Fluor 488



Western blot of human RASSF1A protein in HeLa cell lysate using anti-RASSF1A antibody



Western blot of Human RASSF1A antibody STJA0041537 showing 293T and RASSF1A-transfected lysates