

Anti-SNAP25 antibody (1-206aa) [4E11] (STJA0041539)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-G affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG1k
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	6616
Gene Symbol	SNAP25
UniProt ID	SNP25_HUMAN
Immunogen Region	1-206aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC IF of SNAP25 in U87MG cells stained with STJA0041539 antibody (Alexa Fluor 488 secondary green) and DAPI blue



Western blot of human SNAP25 antibody in mouse brain tissue lysate using anti-mouse HRP detection



Western blot of human SNAP25 protein in 293T cells using anti-SNAP25 antibody