

Anti-GRP94/HSP90B1 antibody (676-803aa) [2H3] (STJA0041542)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF/FACS
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-G affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG2ak
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	7184
Gene Symbol	HSP90B1
UniProt ID	ENPL_HUMAN
Immunogen Region	676-803aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Immunofluorescence staining of human ES cell colony with monoclonal anti-human TRA1 antibody 2H3



Western blot of human GRP94 HSP90B1 protein in ES cell lysate using STJA0041542 antibody



Flow cytometry of Human GRP94 HSP90B1 in U87MG cells with Alexa Fluor 488 detection