

Anti-TPD52L1 antibody (1-131aa) [d1D5] (STJA0041544)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/IHC/FACS
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-G affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG2bk
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	7164
Gene Symbol	TPD52L1
UniProt ID	TPD53_HUMAN
Immunogen Region	1-131aa

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Flow cytometry of TPD52L1 in A431 cells with Alexa Fluor 488 detection



Human gastric cancer tissue immunohistochemistry stained with anti-TPD52L1 antibody using DAB detection