

Anti-ADRM1 antibody (1-407aa) [AT34C2] (STJA0041551)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF/FACS
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-G affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG2bk
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	11047
Gene Symbol	ADRM1
UniProt ID	ADRM1_HUMAN
Immunogen Region	1-407aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC IF of ADRM1 in MCF7 cells stained with STJA0041551 antibody (Alexa Fluor 488 secondary) and DAPI



Western blot of human ADRM1 protein expression across 8 cell lines using anti-human ADRM1 antibody



Flow cytometry of ADRM1 in MCF7 cells with Alexa Fluor 488 staining and IgG isotype control