

Anti-Casein kinase 1 alpha/CSNK1A1 antibody (1-337aa) [AT2E2] (STJA0041567)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-G affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG2bk
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	1452
Gene Symbol	CSNK1A1
UniProt ID	KC1A_HUMAN
Immunogen Region	1-337aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC IF staining of Casein Kinase 1 alpha in HeLa cells with STJA0041567 antibody Alexa Fluor 488 green DAPI blue



Western blot of human Casein Kinase 1 alpha CSNK1A1 in Jurkat cell lysate



Western blot of human Casein Kinase 1 alpha CSNK1A1 in Jurkat cell lysate