

Anti-TBCA antibody (1-108aa) [AT1A5] (STJA0041572)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF/FACS
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-G affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG2bk
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	6902
Gene Symbol	TBCA
UniProt ID	TBCA_HUMAN
Immunogen Region	1-108aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




ICC IF of TBCA in A549 cells stained with STJA0041572 antibody Alexa Fluor 488 secondary DAPI nuclear counterstain




ICC IF of TBCA in HeLa cells stained with STJA0041572 antibody Alexa Fluor 488 secondary DAPI nucleus



Western blot of human TBCA antibody showing LNCaP cell lysate detection with anti-mouse HRP secondary



Western blot of human TBCA antibody on HeLa and K562 cell lysates using anti-mouse HRP detection



Flow cytometry of TBCA protein in HeLa cells using Human TBCA antibody STJA0041572