

Anti-Bcl-2 antibody (1-211aa) [AT1B5] (STJA0041594)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF/FACS
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-G affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG2bk
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	596
Gene Symbol	BCL2
UniProt ID	BCL2_HUMAN
Immunogen Region	1-211aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Immunofluorescence of HeLa cells stained with anti-Bcl-2 antibody and Hoechst 3342



Western blot showing recombinant human Bcl-2 protein detected with anti-human Bcl2 antibody



Flow cytometry of Bcl-2 in HeLa cells using STJA0041594 antibody with Alexa Fluor 488 detection