

Anti-Cytokeratin 18 antibody (79-430aa) [AT2G3] (STJA0041597)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF/IHC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-G affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG2bk
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	3875
Gene Symbol	KRT18
UniProt ID	K1C18_HUMAN
Immunogen Region	79-430aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC IF of Cytokeratin 18 in HeLa cells with Alexa Fluor 488 and DAPI staining



Western blot of human Cytokeratin 18 antibody STJA0041597 in HeLa A549 and MCF7 cell lysates



IHC staining of cytokeatin 18 in paraffin-embedded colorectal cancer tissue using STJA0041597 antibody