

Anti-ARF1 antibody (1-181aa) [AT1B3] (STJA0041609)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF/FACS
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-G affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG1k
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	375
Gene Symbol	ARF1
UniProt ID	ARF1_HUMAN
Immunogen Region	1-181aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




ICC IF of ARF1 in MCF7 cells with STJA0041609 antibody Alexa Fluor 488 and DAPI



ICC IF of ARF1 in U87MG cells stained with STJA0041609 antibody Alexa Fluor 488 secondary and DAPI



Western blot of human ARF1 antibody detecting protein in HeLa HepG2 MCF7 cell lysates and brain tissue



Flow cytometry of human ARF1 in MCF7 cells with Alexa Fluor 488 detection