

Anti-Kallikrein 3/PSA antibody (25-261aa) [AT1D10] (STJA0041628)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF/FACS
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-G affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG1k
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	354
Gene Symbol	KLK3
UniProt ID	KLK3_HUMAN
Immunogen Region	25-261aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC IF of KLK3 in LNCaP cells stained with STJA0041628 antibody Alexa Fluor 488 and DAPI



Western blot of human Kallikrein 3 PSA using anti-KLK3 antibody in PC-3 cell extracts



Flow cytometry of KLK3 in LNCaP cells with Human Kallikrein 3 PSA antibody staining