

Anti-MEMO1 antibody (1-297aa) [AT1E9] (STJA0041630)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/FACS
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-G affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG1k
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51072
Gene Symbol	MEMO1
UniProt ID	MEMO1_HUMAN
Immunogen Region	1-297aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot of human MEMO1 protein in Jurkat cell lysate using anti-MEMO1 antibody



Western blot of human MEMO1 antibody showing U87MG and mouse brain lysate detection



Flow cytometry of MEMO1 in U87MG cells with goat anti-mouse IgG Alexa Fluor 488