

Anti-CNDP2/CPGL antibody (1-475aa) [AT15E5] (STJA0041636)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-G affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG1k
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	55748
Gene Symbol	CNDP2
UniProt ID	CNDP2_HUMAN
Immunogen Region	1-475aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC IF of CNDP2 CPGL in HeLa cells stained with STJA0041636 antibody and Alexa Fluor 488



Western blot of mouse kidney lysate probed with anti-human CNDP2 CPGL antibody



Western blot of human CNDP2 CPGL antibody STJA0041636 showing protein detection in mouse liver and HepG2 lysates