

Anti-SAR1B antibody (1-198aa) [AT1C7] (STJA0041641)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF/FACS
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-A affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG1k
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	51128
Gene Symbol	SAR1B
UniProt ID	SAR1B_HUMAN
Immunogen Region	1-198aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC IF of SAR1B in Hep3B cells with Alexa Fluor 488 and DAPI staining



Western blot showing Human SAR1B antibody detection in mouse muscle tissue lysate using anti-mouse HRP secondary



Flow cytometry of SAR1B in Hep3B cells using human SAR1B antibody STJA0041641