

Anti-MyD88 antibody (1-309aa) [AT22F11] (STJA0041647)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF/FACS
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-A affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG2bk
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	4615
Gene Symbol	MYD88
UniProt ID	MYD88_HUMAN
Immunogen Region	1-309aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC IF of MYD88 in A549 cells stained with STJA0041647 antibody Alexa Fluor 488 secondary green DAPI nucleus blue



Western blot of human MyD88 protein in Raji cell lysate using anti-human MyD88 antibody



Flow cytometry of MYD88 in A549 cells with goat anti-mouse IgG Alexa Fluor 488