

## Anti-ADI1 antibody (1-179aa) [AT27E8] (STJA0041648)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-A affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG1k
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION


Gene ID	<a href="#">55256</a>
Gene Symbol	<a href="#">ADI1</a>
UniProt ID	<a href="#">MTND_HUMAN</a>
Immunogen Region	1-179aa

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC IF of ADI1 in Hep3B cells with STJA0041648 antibody and Alexa Fluor 488 staining



Western blot of human ADI1 antibody in HepG2 cell lysates using anti-mouse HRP detection