

Anti-Hexosaminidase A/HEXA antibody (89-529aa) [AT20F1] (STJA0041656)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/FACS
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-A affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG2aL
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	3073
Gene Symbol	HEXA
UniProt ID	HEXA_HUMAN
Immunogen Region	89-529aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot of human Hexosaminidase A using anti-HEXA antibody on MCF7 cell lysate



Western blot of human HEXA antibody detection in 293T and transfected cell lysates



Flow cytometry of HEXA protein in A549 cells using Human Hexosaminidase A antibody with Alexa Fluor 488 detection