

Anti-CNBP antibody (1-170aa) [AT38F10] (STJA0041658)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/FACS/ICC/IF
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-A affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG2bk
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	7555
Gene Symbol	CNBP
UniProt ID	CNBP_HUMAN
Immunogen Region	1-170aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




ICC IF of CNBP in HeLa cells with DAPI nuclear staining and Alexa Fluor 488 antibody detection




Western blot of Human CNBP antibody STJA0041658 showing protein detection




ICC IF of CNBP in MCF7 cells with DAPI nuclear staining and Alexa Fluor 488 detection



Flow cytometry of Human CNBP antibody STJA0041658 showing fluorescence distribution



Western blot of human CNBP protein in 293T cell lysate detected with CNBP antibody



Flow cytometry of CNBP in MCF7 cells using Human CNBP antibody with Alexa Fluor 488 detection