

Anti-REXO2 antibody (26-237aa) [AT59D9] (STJA0041664)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-A affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG2bk
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	25996
Gene Symbol	REXO2
UniProt ID	ORN_HUMAN
Immunogen Region	26-237aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




ICC IF of REXO2 in 293T cells with DAPI and Alexa Fluor 488 staining




Western blot of Human REXO2 antibody STJA0041664 showing protein detection



ICC IF of REXO2 in HeLa cells with DAPI and Alexa Fluor 488 staining



Western blot of Human REXO2 antibody STJA0041664 showing protein detection



Western blot of human REXO2 protein expression across SW480 HeLa 293T Jurkat K562 and NIH3T3 cell lysates