

Anti-PPM1F antibody (1-454aa) [AT6D11] (STJA0041671)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-A affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG1k
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	9647
Gene Symbol	PPM1F
UniProt ID	PPM1F_HUMAN
Immunogen Region	1-454aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC IF of PPM1F in HeLa cells stained with STJA0041671 antibody Alexa Fluor 488 secondary and DAPI nuclear stain



Western blot of human PPM1F protein in Jurkat cell lysate using anti-PPM1F antibody



Western blot of human PPM1F antibody showing 293T and PPM1F-transfected 293T cell lysates