

Anti-PCBP1 antibody (1-163aa) [AT2A10] (STJA0041676)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-A affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG2ak
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	5093
Gene Symbol	PCBP1
UniProt ID	PCBP1_HUMAN
Immunogen Region	1-163aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC IF of PCBP1 in HeLa cells with Alexa Fluor 488 and DAPI staining



ICC IF of PCBP1 in MCF-7 cells with STJA0041676 antibody and Alexa Fluor 488 detection



Western blot of human PCBP1 protein expression in multiple cell lines using anti-human PCBP1 antibody