

Anti-DNAL1 antibody (1-190aa) [AT29E4] (STJA0041678)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF/FACS
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-A affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG2bk
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	83544
Gene Symbol	DNAL1
UniProt ID	DNAL1_HUMAN
Immunogen Region	1-190aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




ICC IF of DNAL1 in HeLa cells stained with STJA0041678 antibody Alexa Fluor 488 secondary DAPI nucleus



ICC IF of DNAL1 in A549 cells with Alexa Fluor 488 and DAPI staining



Western blot of human DNAL1 antibody in HeLa HepG2 and MCF7 cell lysates using anti-mouse HRP detection



Flow cytometry of DNAL1 in HeLa cells with goat anti-mouse IgG Alexa Fluor 488