

## Anti-RANBP3 antibody (235-445aa) [AT12E11] (STJA0041692)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | ELISA/WB           |
| Host / Source | Mouse              |
| Reactivity    | Human              |

### PRODUCT PROPERTIES


|                     |                                                                                                                                                     |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                                                                          |
| Concentration       | 1 mg/mL                                                                                                                                             |
| Conjugation         | Unconjugated                                                                                                                                        |
| Purification        | By protein-A affinity chromatography                                                                                                                |
| Formulation         | Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol                                                                  |
| Isotype             | IgG1k                                                                                                                                               |
| Storage Instruction | For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles. |

### TARGET INFORMATION

|                  |                             |
|------------------|-----------------------------|
| Gene ID          | <a href="#">8498</a>        |
| Gene Symbol      | <a href="#">RANBP3</a>      |
| UniProt ID       | <a href="#">RANB3_HUMAN</a> |
| Immunogen Region | 235-445aa                   |

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot of human RANBP3 protein in Jurkat cell lysate detected with RANBP3 antibody