

Anti-CACYBP antibody (1-185aa) [AT3G8] (STJA0041697)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-A affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG2ak
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	27101
Gene Symbol	CACYBP
UniProt ID	CYBP_HUMAN
Immunogen Region	1-185aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC IF of CACYBP in PC3 cells with DAPI and anti-CACYBP antibody staining



Western blot of human CACYBP antibody STJA0041697 protein detection



Western blot of CACYBP protein in HeLa 293T Jurkat and A549 cell lysates using anti-human CACYBP antibody