

Anti-CRADD antibody (1-199aa) [AT14G8] (STJA0041707)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-A affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG1ak
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	8738
Gene Symbol	CRADD
UniProt ID	CRADD_HUMAN
Immunogen Region	1-199aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




ICC IF of CRADD in HeLa cells with STJA0041707 antibody and Alexa Fluor 488 detection



ICC IF of CRADD in PC3 cells stained with human CRADD antibody STJA0041707 Alexa Fluor 488 secondary antibody and DAPI nuclear stain



Western blot of human CRADD antibody in A549 HeLa MCF-7 and 293T cell lysates



Western blot of recombinant human CRADD protein detected with anti-human CRADD antibody