

Anti-BMP-8b antibody (264-402aa) [AT13E6] (STJA0041711)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-A affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG2ak
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	656
Gene Symbol	BMP8B
UniProt ID	BMP8B_HUMAN
Immunogen Region	264-402aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC IF of BMP-8b protein in HaLe cells stained with STJA0041711 antibody Alexa Fluor 488 secondary DAPI nuclear counterstain



Western blot of human BMP-8b antibody in 293T cells using STJA0041711



Western blot of recombinant human BMP-8b protein using anti-BMP-8b antibody at 1:1000 dilution