

Anti-CD14 antibody (20-349aa) [AT87H7] (STJA0041718)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/FACS
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-A affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG2bk
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	929
Gene Symbol	CD14
UniProt ID	CD14_HUMAN
Immunogen Region	20-349aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot of human CD14 antibody STJA0041718 in HeLa and A549 cell lysates



Flow cytometry of CD14 in THP-1 cells stained with STJA0041718 antibody



Flow cytometry of CD14 in PBMC cells using STJA0041718 antibody with Alexa Fluor 488 detection