

Anti-Aldolase A/ALDOA antibody (1-364aa) [AT3F9] (STJA0041720)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF/FACS
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-A affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG2ak
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	226
Gene Symbol	ALDOA
UniProt ID	ALDOA_HUMAN
Immunogen Region	1-364aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC IF of ALDOA in A549 cells with STJA0041720 antibody and Alexa Fluor 488 detection



Western blot showing Human ALDOA antibody detection in recombinant protein and cell lysates



Flow cytometry of ALDOA in A549 cells with Alexa Fluor 488 staining